

Labware

Click on your product of interest for more information.

Sealing Film Guide




Product	SealPlate® Films	SealPlate® Ministrips™ Films	EZ-Pierce™ Films	Zone-Free™ Films	X-Pierce™ Films	AbsorbMax™ Films	ThermaSeal™ Films	ThermaSeal RT™ Films	AlumaSeal II™ Films	AlumaSeal 96™ Films	AlumaSeal 384™ Films	AeraSeal™ Films
Application	ELISA, Incubation, Storage	ELISA, Incubation, Storage	Processing, Autosampling	Processing, Autosampling	Processing, Autosampling	Fluorescence, Storage	PCR, Storage	Real-time PCR	PCR, Cold Storage	PCR, Storage	PCR, Storage	Cell Culture, Tissue Culture
Special Properties	Eliminates edge effects	2x8-well strips	Easily pierceable, chemical resistant, conforms to plate	Easily pierceable, chemical resistant, adhesive-free area over each well	Pre-cut, pierceable, resealing, printed well guides	Black, light absorbing	Chemical resistant, good moisture barrier	Best optical clarity	Chemical resistant, light blocking, good moisture barrier	Chemical resistant, light blocking, good moisture barrier	Chemical resistant, light blocking, good moisture barrier	Breathable
Minimum Temp.	-40 °C	-40 °C	-40 °C	-40 °C	-40 °C	-40 °C	-40 °C	-40 °C	-80 °C	-80 °C	-80 °C	-40 °C
Maximum Temp.	+120 °C	+120 °C	+90 °C	+90 °C	+90 °C	+80 °C	+125 °C	+120 °C	+120 °C	+150 °C	+150 °C	+80 °C
Available Sterile	Yes	Yes	Yes	Yes	Yes	Yes	Yes		Yes			Yes
Pierceable			Yes	Yes	Yes				Yes	Yes	Yes	Yes
RNase/DNase Free							Yes	Yes	Yes	Yes	Yes	
Material	Polyester	Polyester	Polyethylene	Polyethylene/Polypropylene	Vinyl	Vinyl	Polyethylene	Polyester	Aluminum	Aluminum	Aluminum	Rayon
NonSterile	Z369659	Z707384	Z721581	Z721646	Z722502	Z722537	Z369675	2 mil thickness Z707465 5 mil thickness Z707473	A2350	Z721549	Z721557	Z707430
Sterile	Z369667	Z707392	Z721603	Z721654	Z722529	Z722545	Z369683		Z707422			A9224

New!




Available with ColorTab™ End Tabs

- Z721417
- Z721425
- Z721433
- Z721441
- Z721468
- Z721476
- Z721484




Available printed well guides and ColorTab™ End Tabs

- Z721611
- Z721638




New!



Also available in MiniStrip™ Format Non-Sterile Only
Z721492

ThermaSeal RT2RR™ for Raised-Rim Plates Non-Sterile Only
Z722553



Equivalent to CLS6570

Equivalent to CLS6569

