

HandyStep® S

Don't miss our special price for a repetitive pipette HandyStep® S valid from 1st January until 31st March 2021!

HandyStep® S – the repetitive pipette for so many applications!

Ideal for routine or research applications in diagnostics, molecular biology, environmental analysis, and many more. HandyStep® S and PD-Tips work on the direct displacement principle. This offers the highest-precision dispensing of liquid media with high viscosity, high density, or high vapor pressure. Direct displacement permits contamination-free operation, since no aerosols are formed.

Your benefits:

- Serial pipetting made easy
- Works with BRAND® PD-Tips and third-party tips
- Robust, stable and durable



Save 25% off list price

Description	Cat. No.
Repetitive pipette HandyStep® S, DE-M marking, shelf/rack mount, 3 PD-Tips: 0.1 ml, 1 ml and 10 ml	Z692867-1EA

Please use promo code **UFN** to receive a 25% discount!

Order must reach us between 1st January and 31st March 2021. Valid in Western Europe only (Austria, Germany, Switzerland & Liechtenstein, United Kingdom, Ireland, Norway, Sweden, Finland, Denmark, France, BeNeLux, Spain, Italy, Portugal)

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

© 2020 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck and the vibrant M are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

Lit Code MK_FL7271EN 34344 12/2020

