



# Air Gap

## General Description

The air gap is a small plastic device that is placed on tubing going from the water system to the drain. This creates an unobstructed vertical space, the “air gap” which will prevent contaminants in the drain from going back into the system.

## Key Benefits

- Provides cost-effective protection of the water system and its accessories.
- Prevents bacteriological contamination from the drain.

## Ordering Information

Description	Cat. No.
Air Gap Device	AIRGAP001



## To place an order or receive technical assistance

In the U.S. and Canada, call toll-free  
1-800-645-5476

For Technical Service, please visit:  
[emdmillipore.com/techservice](http://emdmillipore.com/techservice)

For other countries across Europe  
and the world, please visit:  
[emdmillipore.com/offices](http://emdmillipore.com/offices)

[milliporesigma.com](http://milliporesigma.com)

The life science business of Merck KGaA, Darmstadt, Germany operates as  
MilliporeSigma in the U.S. and Canada.

© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. MilliporeSigma, Milli-Q, and the vibrant M are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MS\_DS1252EN  
2017 - 08107  
1/2018