

Product Information

Mineral Oil, PCR Reagent

Catalog Number **M8662**

Product Profile

DNase, RNase and protease: None detected
Suitable for use in the Polymerase Chain Reaction (PCR).

PCR Suitability

100 µl of mineral oil was overlaid on a 100 µl PCR reaction containing: 10 mM Trizma®-HCl, pH 8.3 at 25 °C, 50 mM KCl, 1.5 mM MgCl₂, 0.001% (w/v) gelatin, each dNTP at 200 µM, primers defining an approximately 500 base pair region of λ DNA at 1.0 :M each, λ DNA template at 1 ng/100 µl, and Taq DNA Polymerase at 2.5 units/100 µl. The reaction underwent 25 cycles of 94 °C to denature the double stranded DNA, 55 °C to anneal the DNA segments, and 72 °C to extend the DNA segments. A single band of approximately 500 base pairs was visualized following electrophoresis of the reaction product in a 1.5% agarose gel.

Endonuclease-Exonuclease

One µg of λ Hind III fragments was incubated for 16 hours at 37 °C with 25 µl mineral oil overlay on a 50 µl reaction mixture containing 30 mM Trizma®-HCl, pH 7.8, 50 mM NaCl and 10 mM MgCl₂. No degradation of the DNA fragments was detected following agarose gel electrophoresis. Detection limit: Degradation of 10% of the DNA substrate is detectable.

Endonuclease (Nickase)

One µg of pBR322 DNA was incubated for 16 hours at 37 °C with 25 µl mineral oil overlay on a 50 µl reaction mixture containing 30 mM Trizma®-HCl, pH 7.8, 50 mM NaCl and 10 mM MgCl₂. No conversion of the covalently closed circular DNA to the nicked or linear form was observed following agarose gel electrophoresis. Detection limit: Conversion of 1% of the DNA substrate is detectable.

RNase

Two µg of transfer RNA were incubated for 16 hours at 37 °C with 25 µl mineral oil overlay on a 50 µl reaction mixture containing 30 mM Trizma®-HCl, pH 7.8, 50 mM NaCl and 10 mM MgCl₂. No degradation of the tRNA was detected following polyacrylamide gel electrophoresis. Detection limit: Degradation of 10% of the tRNA substrate is detectable.

Protease

50 µl of 0.5% FITC-casein reaction mixture was overlaid with 10 µl of mineral oil and incubated for 1 hour at 37 °C. Liberated FITC equivalents are quantitated fluorometrically. Detection limit: 8.5 X 10⁻⁶ µmoles of FITC released per minute.

Modified from: Twining, S.S., Analytical Biochemistry **143**, 30-34 (1984)

Information about licenses to PCR can be obtained from The Perkin-Elmer Corporation or Roche Molecular Systems, Inc.

PHC 11/10-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.

Please see reverse side of the invoice or packing slip.