



**QUALITY CONTROL RANGES**  
**MILLIPLEX® MAP Human Fibrosis Panel 1**  
**Catalog # HFIB1-100K**  
**Control Catalog # HFIB1-6100**

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
IL-18BPα	Control I	HFIB1-102	66 - 137	pg/mL
	Control II	HFIB1-202	334 - 694	pg/mL
COMP	Control I	HFIB1-102	251 - 522	pg/mL
	Control II	HFIB1-202	1336 - 2775	pg/mL
Uteroglobin	Control I	HFIB1-102	139 - 289	ng/mL
	Control II	HFIB1-202	729 - 1514	ng/mL

**Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit.  
Quality Control values in culture media are not tested.**

HFIB1-6100: HFIB1-102.202 (13 January 2023)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) [emdmillipore.com](http://emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA