



204862




16-Square Well ProPlate Slide Module Footprint

Number of wells	Item No.	Quantity Per Pack
16	204862	2




204860




16-Round Well ProPlate Slide Module Footprint

Number of wells	Item No.	Quantity Per Pack
16	204860	2




246865




64-Well ProPlate Slide Module Footprint

Number of wells	Item No.	Quantity Per Pack
64	246865	2



246861


1-Well ProPlate Slide Module Footprint



Number of wells	Item No.	Quantity Per Pack
1	246861	2

246862


2-Well ProPlate Slide Module Footprint



Number of wells	Item No.	Quantity Per Pack
2	246862	2


246863

3-Well ProPlate Slide Module Footprint




Number of wells	Item No.	Quantity Per Pack
3	246863	2

246864




4-Well ProPlate Slide Module Footprint

Number of wells	Item No.	Quantity Per Pack
4	246864	2




246868



8-Well ProPlate Slide Module Footprint

Number of wells	Item No.	Quantity Per Pack
8	246868	2




204838

Stainless Steel Spring Clips




Spring Clip  
End View

Stainless steel spring clips are recommended for use with glass slides which have un-polished edges.



Slide module assembled  
using Spring Clips  
End View



Slide module assembled  
using Spring Clips  
Top View

Item  
No.


Quantity  
Per  
Pack

204838

10

204874

Seal Strips and Applicator



Item  
No.

Quantity  
Applicators  
Per Pack


Quantity  
Seal Strips  
Per Pack

204874

1

50

ProPlate Tray Set



204860 Shown

Item  
No.

Number  
of wells

Quantity  
Per  
Pack

204860

16 SQ

1

247880

16 RND

1

247881

1

1

247882

2

1

247883

3

1

247884

4

1

247888

8

1

247885

64

1

Each ProPlate Tray Set includes: 4 Slide modules, 8 Snap clips, 1 Tray and Cover, 10 seal strips and one applicator.

Tray conforms to SBS standards for footprint (81.7mm width x 123.7mm length) and well spacing.