

don't use printout for evaluation!!!





1.14410
MQuant®
Silicate Test
Silicat-Test
Test Silicates
Test Silicatos

Si

Color card
 Farbkarte
 Carte colorimétrique
 Tarjeta colorimétrica



© 2021



Merck KGaA, 64271 Darmstadt, Germany,
 Tel. +49(0)6151 72-2440
www.sigmaldrich.com/mobile-laboratory

EMD Millipore Corporation, 400 Summit Drive
 Burlington MA 01803, USA, Tel. +1-978-715-4321

Sigma-Aldrich Canada Co. or Millipore (Canada) Ltd.
 2149 Winston Park, Dr. Oakville, Ontario, L6H 6J8
 Phone: +1 800-565-1400

7.69144.1007-XXXXXX 09/2021

0 **0.01** **0.02** **0.04** **0.06** **0.08** **0.10** **0.15** **0.20** **0.25** mg/l Si



20 ml



3 x Si-1



pH 1.2 - 1.6



3 min




3 x Si-2



10 x Si-3



2 min



0 **0.021** **0.043** **0.086** **0.13** **0.17** **0.21** **0.32** **0.43** **0.53** mg/l SiO₂

