

shelf life information data sheet

Product Family: MILLICELL® HANGING CELL CULTURE INSERT

shelf life information

Shelf Life: 36 Months

Expiration Date: See product label

Storage Conditions: Ambient temperature

Applicable Catalog Number(s): MCHT06H48, MCRP06H48, MCSP06H48, MCMP06H48, MCEP06H48,
MCHT12H48, MCRP12H48, MCSP12H48, MCMP12H48, MCEP12H48,
MCHT24H48, MCRP24H48, MCSP24H48, MCMP24H48, MCEP24H48

quality attributes studied

Visual Inspection

Status of Membrane seal area

Membrane Pore Size

Membrane Thickness

Dose of Gamma Irradiation

The stated expiration date designates the time prior to which a lot of the product is expected to remain within the approved shelf life specification if stored under defined conditions, and after which must not be put into use.

Not all products listed have undergone stability testing to determine the products' shelf lives. The design of the stability schedule is such that only select product(s) within a product family are tested. The design assumes that the products selected are representative of the family as they utilize the same materials of construction, production processes, and packaging components as the products within the family.

Updates are made to this data sheet on an as-needed basis; the date of signature on this document reflects the most recent revision.



28-Oct-19

Thomas Brackett
Supplier Quality Engineering Manager
Third Party Finished Goods Quality

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

