

Product Information

M2 Medium with HEPES, without penicillin and streptomycin

Catalog Number **M7167**
Storage Temperature 2–8 °C

Product Description

M2 and M16 Medium are common media for *in vitro* culture of preimplantation stage embryos. It is a modified Krebs-Ringer bicarbonate solution, which is very similar to Whitten's Medium. M16 contains pyruvate and lactate as energy sources since preimplantation embryos do not utilize glucose efficiently. M2 Medium is a further modification of M16 that substitutes HEPES buffer in place of some of the bicarbonate. M2 is used for collecting and handling embryos for prolonged periods outside a CO₂ incubator.

Biological Performance: This product is tested for its ability to support the development of one-cell mouse embryos to expanded blastocysts. Minimum requirement is 80% development to blastocyst.

Components	g/L
Calcium Chloride · 2H ₂ O	0.25137
Magnesium Sulfate (anhydrous)	0.1649
Potassium Chloride	0.35635
Potassium Phosphate, Monobasic	0.162
Sodium Bicarbonate	0.35
Sodium Chloride	5.53193
Albumin, Bovine Fraction V	4.0
D-Glucose	1.0
HEPES · Na	5.42726
Phenol Red · Na	0.0106
Pyruvic Acid · Na	0.0363
DL-Lactic Acid · Na	2.95

Supplement with 0.06 g/L potassium penicillin-G and 0.05 g/L streptomycin sulfate.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Store the liquid medium at 2–8 °C in the dark. The nature of supplements added may affect storage conditions and shelf life of the medium. Product label bears expiration date. Deterioration of the liquid medium may be recognized by any of the following:

1. pH change
2. cloudy appearance
3. color change.

Note: Slight precipitate or particulate matter may be seen throughout the solution, which does not affect the performance of the product.

References

1. Whitten, W.K., *Adv. Biosci.*, **6**, 129-141 (1971).
2. Quinn, P. *et al.*, *J. of Reprod. Fertil.*, **66**, 161-168 (1982).
3. Hogan, B. *et al.*, *Manipulating the Mouse Embryo*. Cold Spring Harbor Laboratory. Cold Spring Harbor, NY. Section G, pp. 249-252 (1986).

GS,LM,MAM 02/10-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.