






## Product Information

### SIGMA QUALITY CONTROL TEST PROCEDURE

#### Enzymatic Assay of ALDOLASE<sup>1</sup> (EC 4.1.2.13)

##### PRINCIPLE:



Abbreviations used:

G 3-P = Glyceraldehyde 3-Phosphate

DHAP = Dihydroxyacetone Phosphate

TPI = Triosephosphate Isomerase

$\alpha$ -GDH =  $\alpha$ -Glycerophosphate Dehydrogenase

$\beta$ -NADH =  $\beta$ -Nicotinamide Adenine Dinucleotide, Reduced Form

$\beta$ -NAD =  $\beta$ -Nicotinamide Adenine Dinucleotide, Oxidized Form

**CONDITIONS:** T = 25°C, pH = 7.4, A<sub>340nm</sub>, Light path = 1 cm

**METHOD:** Continuous Spectrophotometric Rate Determination

##### REAGENTS:

- A. 100 mM Tris HCl Buffer, pH 7.4 at 25°C.  
(Prepare 250 ml in deionized water using Trizma Base, Sigma Prod. No. T-1503. Adjust to pH 7.4 at 25°C with 1 M HCl.)
- B. 58 mM Fructose 1,6-Diphosphate Solution (F 1,6-DP)  
(Prepare 1 ml in deionized water using D-Fructose 1,6-Diphosphate, Tetra(cyclohexylammonium) Salt, Sigma Prod. No. F-0752.)
- C. 4.0 mM  $\beta$ -Nicotinamide Adenine Dinucleotide, Reduced Form Solution ( $\beta$ -NADH)  
(Prepare 2 ml in cold deionized water using  $\beta$ -Nicotinamide Adenine Dinucleotide, Reduced Form, Disodium Salt, Sigma Prod. No. N-8129 or dissolve the contents of one 5 mg vial of  $\beta$ -Nicotinamide Adenine Dinucleotide, Reduced Form, Disodium Salt, Sigma Stock. No. 340-105, in the appropriate volume of deionized water.)
- D.  $\alpha$ -Glycerophosphate Dehydrogenase/Triosephosphate Isomerase Enzyme Solution ( $\alpha$ -GDH/TPI)  
(Immediately before use, prepare a solution containing 50  $\alpha$ -GDH units/ml of  $\alpha$ -Glycerophosphate Dehydrogenase/Triosephosphate Isomerase, Sigma Prod. No. G-6755, in cold deionized water.)

#### Enzymatic Assay of ALDOLASE<sup>1</sup> (EC 4.1.2.13)

**REAGENTS:** (continued)

- E. Aldolase Enzyme Solution  
(Immediately before use, prepare a solution containing 0.25 - 0.50 unit/ml of Aldolase in cold Reagent A.)

**PROCEDURE:**

Pipette (in milliliters) the following reagents into suitable cuvettes:

	<u>Test</u>	<u>Blank</u>
Reagent A (Buffer)	2.60	2.60
Reagent B (F 1,6-DP)	0.10	0.10
Reagent C ( $\beta$ -NADH)	0.10	0.10
Reagent D ( $\alpha$ -GDH/TPI)	0.10	0.10

Mix by inversion and equilibrate to 25°C. Monitor the  $A_{340\text{nm}}$  until constant, using a suitably thermostatted spectrophotometer. Then add:

Reagent A (Buffer)	-----	0.10
Reagent E (Enzyme Solution)	0.10	-----

Immediately mix by inversion and record the decrease in  $A_{340\text{nm}}$  for approximately 5 minutes. Obtain the  $\Delta A_{340\text{nm}}/\text{min}$  using the maximum linear rate for both the Test and Blank.

**CALCULATIONS:**

$$\text{Units/ml enzyme} = \frac{(\Delta A_{340\text{nm}}/\text{min Test} - \Delta A_{340\text{nm}}/\text{min Blank})(3)(\text{df})}{(2)(6.22)(0.1)}$$

3 = Total volume (in milliliters) of assay

df = Dilution factor

2 = 2 moles of  $\beta$ -NADH converted to 2 moles of  $\beta$ -NAD per mole of Fructose 1,6-Diphosphate

6.22 = Millimolar extinction coefficient of  $\beta$ -NADH at 340 nm

0.1 = Volume (in milliliter) of enzyme used

**Enzymatic Assay of ALDOLASE<sup>1</sup>**  
**(EC 4.1.2.13)**

**CALCULATIONS:** (continued)

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

**UNIT DEFINITION:**

One unit will convert 1.0  $\mu$ mole of fructose 1,6-diphosphate to dihydroxyacetone phosphate and glyceraldehyde 3-phosphate per minute at pH 7.4 at 25°C.

**FINAL ASSAY CONCENTRATION:**

In a 3.00 ml reaction mix, the final concentrations are 90 mM Tris, 1.9 mM fructose 1,6-diphosphate, 0.13 mM  $\beta$ -nicotinamide adenine dinucleotide, 5 units  $\alpha$ -glycerophosphate dehydrogenase/triosephosphate isomerase (based on  $\alpha$ -glycerophosphate dehydrogenase units) and 0.025 - 0.050 unit aldolase.

**REFERENCE:**

Bergmeyer, H.U. (1974) *Methods of Enzymatic Analysis*, Second Edition, Volume I, 430

**NOTES:**

1. This enzyme assay is not to be used to assay Aldolase, from *Staphylococcus aureus*, Sigma Prod. No. A-2548, Aldolase, insoluble enzyme attached to polyacrylamide from Rabbit Muscle, Sigma Prod. No. A-1386, and Aldolase from Baker's Yeast, Sigma Prod. No. A-9562.
2.  $\alpha$ -Glycerophosphate Dehydrogenase Unit Definition: One unit will convert 1.0  $\mu$ mole of dihydroxyacetone phosphate to  $\alpha$ -glycerophosphate per minute at pH 7.4 at 25°C.
3. Triosephosphate Isomerase Unit Definition: One unit will convert 1.0  $\mu$ mole of D-glyceraldehyde 3-phosphate to dihydroxyacetone phosphate per minute at pH 7.6 at 25°C.
4. All product and stock numbers, unless otherwise indicated, are Sigma product and stock numbers.

**Sigma warrants that the above procedure information is currently utilized at Sigma and that Sigma products conform to the information in Sigma publications. Purchaser must determine the suitability of the information and products for its particular use. Upon purchase of Sigma products, see reverse side of invoice or packing slip for additional terms and conditions of sale.**