



## GOAT ANTI MOUSE IgG AFFINITY PURIFIED SECONDARY ANTIBODY

<b>CATALOG NUMBER:</b>	AP124
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	2 mg
<b>CONCENTRATION:</b>	XX mg/mL
<b>SPECIFICITY:</b>	Reacts with Mouse IgG gamma chain, as well as the light chains from all mouse immunoglobulin classes. Monospecific on IEP against whole mouse serum.
<b>IMMUNOGEN:</b>	Highly purified Mouse IgG isolated from pooled normal mouse serum.
<b>PURITY:</b>	Single band on IEP at a concentration of 10 mg/mL when run against rabbit anti goat whole serum.
<b>FORMAT:</b>	Liquid in 0.02M Phosphate Buffer, 0.25M Sodium Chloride and 0.1% Sodium Azide. Filtered through a 0.45 micron membrane filter.
<b>STORAGE:</b>	Store at +2-8°C in undiluted aliquots for up to six months.

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200  $\mu$ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC  
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.



**upstate** | CHEMICON | *Linco*

THE EXPERTISE OF UPSTATE®, CHEMICON® AND LINCO®  
IS NOW A PART OF MILLIPORE