

## Technical Bulletin

# Albumin Quick Test Strips

**Catalog Number MAS009**

## Product Description

Albumin is the most abundant plasma protein in humans. Albumin plays many important physiological roles, including maintenance of colloid osmotic pressure and binding of key substances such as long-chain fatty acids, bile acids, bilirubin, hematin, calcium, and magnesium. Serum albumin is a reliable prognostic indicator for liver disease, nephritic syndrome, malnutrition, protein-losing enteropathies, as well as morbidity and mortality. In addition, elevated levels of serum albumin are associated with dehydration.

The Albumin Quick Test Strips are based on the immobilization of an albumin sensitive chromogen onto a pad. Albumin will react with the chromogen and produce a colored product. The intensity of the product color is directly proportional to the albumin concentration in the sample. The semi-quantitative detection range of the kit is 0-0.4 g/dL or 0-60  $\mu$ M (undiluted) albumin.

This kit is suitable for the detection of albumin in serum, plasma, urine, etc.

## Components

The kit is sufficient for 10 tests

- Albumin Test Strips (10 strips) 1 Each  
Catalog Number MAS009A
- Dilution Tubes 10 Each  
Catalog Number MAS009B

## Equipment Required but Not Provided

- Pipetting devices and accessories

## Precautions and Disclaimer

For Research Use Only. Not for use in diagnostic procedures. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

## Storage/Stability

The kit is shipped at room temperature. Store components at 2-8 °C. Keep strips dry and out of direct sunlight.

## Preparation Instructions

### Sample Preparation

Most urine samples can be tested undiluted.

For most serum or plasma samples, it is recommended to dilute samples 100-fold.

For samples requiring dilution, dilute samples into Dilution Tubes according to Table 1. Mix thoroughly.

**Table 1.**  
Sample Dilutions

Dilution	Purified Water	Sample
100-Fold	495 $\mu$ L	5 $\mu$ L

## Procedure

1. Transfer 20  $\mu$ L of diluted or undiluted Sample directly onto the reaction pad of an Albumin Test Strip.
2. Allow color to develop on strip for 30 seconds.
3. Compare the color of the reaction pad of the strip with the provided Albumin Chart. Multiply the concentration in the chart by the Sample dilution factor used (i.e. 100) to determine the Albumin concentration in the original Sample.

---

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

## Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices).

## Technical Service

Visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

## Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).