



QUALITY CONTROL RANGES

MILLIPLIX[®] MAP Human Complement Magnetic Bead Panel 1
Kit Catalog# HCMP1MAG-19K
Control Catalog # HCMP1-6019-1

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
Complement C2	Control 1	HCM1-105	10 – 20	ng/mL
	Control 2	HCM1-205	38 – 79	ng/mL
Complement C4b	Control 1	HCM1-105	9.2 – 19	ng/mL
	Control 2	HCM1-205	37 – 77	ng/mL
Complement C5	Control 1	HCM1-105	19 – 39	ng/mL
	Control 2	HCM1-205	77 – 160	ng/mL
Complement C5a	Control 1	HCM1-105	29 – 61	pg/mL
	Control 2	HCM1-205	118 – 244	pg/mL
Complement Factor D	Control 1	HCM1-105	0.43 – 0.89	ng/mL
	Control 2	HCM1-205	1.8 – 3.7	ng/mL
MBL	Control 1	HCM1-105	0.93 – 1.9	ng/mL
	Control 2	HCM1-205	3.9 – 8.0	ng/mL
Complement Factor I	Control 1	HCM1-105	5.1 – 10	ng/mL
	Control 2	HCM1-205	21 – 44	ng/mL

HCMP1-6019-1: HCM1-105.205 (02 February 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.