


Monoclonal Anti- γ -Tubulin, Clone GTU-88 (T 6557)



Immunoblot

Whole cell extract of chicken fibroblasts was separated on SDS-PAGE and blotted with Monoclonal Anti- γ -Tubulin (Product No. T 6557). The antibody was developed with Goat anti-Mouse IgG, Alkaline Phosphatase conjugate (Product No. A 1293) and NBT/BCIP substrate.

Lane A: Antibody dilution 1:20,000

Lane B: Negative control (only secondary antibody)