

**QUALITY CONTROL RANGES**  
**Mouse Immunoglobulin Isotyping Panel Catalog # MGAMMAG-300K**  
**Control Catalog # 43-008L Lot # MIG-11K**  
**Standard Catalog # 47-300 Lot # 300-11K**



Analyte	Standard Value	Analyte	Standard Value
IgA	5358 ng/mL	IgG2b	2231 ng/mL
IgG1	2233 ng/mL	IgG3	1431 ng/ml
IgG2a	2200 ng/mL	IgM	1305 ng/mL

<u>Analyte</u>	<u>Expected Range</u>	<u>Units</u>
IgA	265 - 551	ng/mL
IgG1	105 - 219	ng/mL
IgG2a	122 - 254	ng/mL
IgG2b	95 - 198	ng/mL
IgG3	59 - 123	ng/mL
IgM	53 - 110	ng/mL

43-008L: MIG-11K (21 February 2020)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • [www.emdmillipore.com](http://www.emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.