



## QUALITY CONTROL RANGES

**MILLIPLEX<sup>®</sup> MAP Human Sepsis Magnetic Bead Panel 1 Kit**  
**Kit Catalog# HSP1MAG-63K**  
**Control Catalog # LHSP-6063**  
**Lot# SPS-115 and SPS-215**

**Standard Catalog # LHSP-8063**  
**Lot# 8063-15K**

Analyte	Standard Value
MIF	74,500 pg/mL

**\*\*Please see reverse side for Quality Control Ranges\*\***

**Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.**

LHSP-6063: SPS-115.215 (16 November 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) [www.emdmillipore.com](http://www.emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.



## QUALITY CONTROL RANGES

**MILLIPLEx<sup>®</sup> MAP Human Sepsis Magnetic Bead Panel 1 Kit**  
**Kit Catalog# HSP1MAG-63K**  
**Control Catalog # LHSP-6063**  
**Lot# SPS-115 and SPS-215**

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
MIF	Control I	490 – 1018	pg/mL
	Control II	4111 – 8539	pg/mL
ICAM-1	Control I	1674 – 3476	pg/mL
	Control II	18039 – 37465	pg/mL
Fas Ligand	Control I	300 – 623	pg/mL
	Control II	3112 – 6464	pg/mL
Soluble Fas	Control I	234 – 486	pg/mL
	Control II	2515 – 5223	pg/mL
VCAM-1	Control I	1441 – 2992	pg/mL
	Control II	14455 – 30022	pg/mL
Total PAI-1	Control I	293 – 608	pg/mL
	Control II	3147 – 6535	pg/mL

**Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.**

LHSP-6063: SPS-115.215 (16 November 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) [www.emdmillipore.com](http://www.emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.