

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Duolink® In Situ Mounting Medium with DAPI

Product Number:

DUO82040

TEST	Specification
------	---------------

Suitability

Suitable for use as a nuclear stain when compared to a control.

Suitable

Specification: PRD.1.ZQ5.1000077758

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.